



Autonomous Medical Officer Support (AMOS) ISS Technology Demonstration: Enabling Earth- Independent Procedure Guidance

Human Research Program
Exploration Medical Capability Element

February 9, 2023

Doug Ebert¹

Vicky Byrne¹

Marlei Walton²

Ashot Sargsyan¹

¹ KBR; ² NASA JSC


“Expanding the Boundaries of Space Medicine and Technology”

- **Problem statement: exploration context**
- **AMOS description**
- **Tech demo**
 - Procedures
 - Results
- **Conclusion**
- **Forward work**






Current



Proposed





- **Fulfills three functions:**
 1. Preflight training
 2. Onboard just-in-time (JIT) or refresher training
 3. Guidance for autonomous medical procedures
- **Intuitive menu-driven, human-centered, modular design**
- **Coded for flexibility and compatibility with ISS systems**
- **Current modules include kidney and urinary bladder ultrasound**
- **Integrated use tracking and evaluation features**
- **Platform technology - extensible to other subject areas**

Web version of AMOS software available at:



<http://comfort2.eastus.cloudapp.azure.com/View/Index/5>

- **Enhance telemedicine capability; introduce asynchrony**
- **Enable autonomous procedure performance**
- **Streamline the process of skill management**
- **Reduce preflight training load**
- **Repeated similar ISS Demos could reduce exploration kidney stone and urinary retention risks (improved monitoring)**

Probe Handling: Reference marker

One side (edge) of the probe has a reference marker (a ridge you can see and feel, and an indentation as shown). This elevated ridge on one side of the ultrasound probe relates to one side of the ultrasound screen. The reference marker orients the ultrasound probe to the ultrasound image display. The respective side of the image (the marker side) is identified by the screen indicator.

Each single ultrasound image shows a 2D cut through tissue under the probe (as shown in right Figure with probe and skin on top, marker side on the left screen indicator).



Representation of a 2-D ultrasound imaging beam




The probe marker is oriented towards the subject's head ('up'). Therefore, the top of the bladder is on the left (indicator) side of the image in this vertical view of a bladder.



Right Kidney: Finding the Kidney and the Starting View

STEP 3 OF 11 OPTIMIZE THE RIGHT KIDNEY STARTING VIEW:

- Apply probe at zone purple 2 (grid in Figure to right) with probe marker towards the subject's head. ([Review probe handling in Foundation section](#))
- Look for the kidney (Note that kidneys move with respiration).
- If not found at once, use rapid meandering up, down, and towards the side of the body. Once the kidney is found, reduce speed and range of probe movement.
- Apply more gel, reapply probe to same location, and adjust the gain as needed.
- Slowly move the probe to hover over the kidney, i.e., bring the kidney to the middle of the image with small adjustments.
- Tune probe position for longest image (left-to-right dimension on the screen) to include upper pole (left on image) and lower pole (right on image) ([Review in Foundation section](#)).
- Translate probe away from spine while keeping the kidney in longest view, optimizing for the widest image (up- down dimension on the screen).
- Practice tilt-sweeping throughout the volume of the kidney. Longest and widest view is the Starting View ([Review in Foundation section](#)).



Ideal Starting View: Note the interruption of the darker tissue at the bottom center of the image, an indication of the longest, widest view.



Finding The Right Kidney and Its Longest, Widest Section. Press play to start video.

PROCEDURE MILESTONE: You have arrived at the Starting View for right kidney imaging- ready to begin right kidney imaging and data collection.

- Flown as ExMC/HRP payload – supported by ROI
- AMOS installed with queso service pack (3/31 - 4/2/2020)
- Demos occurred 4/9/2020 (observed remotely - COVID-19) & 6/29/2022
- Crew Procedures:
 - BEFORE scanning, operator reviews AMOS “Foundations” material
 - Deploy Ultrasound2 per nominal procedures
 - Set up cabin video
 - Open AMOS application
 - Perform bladder scan (full), perform bladder scan (post void)
 - Perform kidney scan (right and left)

- ✓ **Installation / use of AMOS in the ISS server environment (IT objective)**
- ✓ **Use of AMOS for autonomous imaging in operational setting (telemedicine proof of concept objective)**
 - ✓ Obtain click tracking data
 - ✓ Collect ultrasound images - urinary bladder and kidneys (irrespective of quality or compliance)
- ✓ **Collect crew feedback – AMOS implementation (telemedicine and AMOS improvement objective)**
- ✓ **Record in-cabin and ultrasound scanhead video data – technical and human factors assessment (telemedicine and AMOS improvement objective)**
- ✓ **Collect crew feedback – integrated training and procedure support (training, telemedicine, and AMOS improvement objective)**

	Demo 1	12	PROCEEDURE EXECUTION (0-3)	INCLUSION OF TARGET ORGAN VOLUME (0-3)	IMAGING PLANE ORIENTATION (0-3)	SWEEP SPEED ADEQUACY (0-3)	AVERAGE IMAGE QUALITY (0-3)	MEASURABILITY (0-3)	EFFECTIVE IMAGING CONTENT (0-3)	SUCCESS SCORE (0-3)	Y	SUCCESS (Y/N)
	Demo 2	12*	2.5	3.0	3.0	3.0	3.0	2.3	2.5	2.4	Y	

* Three instances were not included in the analysis; operator self-corrected labeling and replaced these images

Demo 2 – Bladder Results (N=1 operator/subject)

	NUMBER OF INSTANCES (CINE-LOOPS)	PROCEDURE EXECUTION (0-3)	INCLUSION OF TARGET ORGAN VOLUME (0-3)	IMAGING PLANE ORIENTATION (0-3)	SWEEP SPEED ADEQUACY (0-3)	AVERAGE IMAGE QUALITY (0-3)	MEASURABILITY (0-3)	EFFECTIVE IMAGING CONTENT (0-3)	SUCCESS SCORE (0-3)	SUCCESS (Y/N)
	3	3	3	3	3	3	3	3	3	Y
	3	3	3	3	3	3	3	3	3	Y
	3	2	3	3	3	1	1	2	1.5	N
	3	2	3	3	3	2	2	2	2	Y
TOTAL NUMBERS AND AVERAGE SCORES		12	2.5	3.0	3.0	2.3	2.3	2.5	2.4	Y

- **Operator initiated “end exam” as written in the procedure**
 - Instructions to skip “end exam” were in the execute note (different location)
 - Operator feedback: include an “if performing an additional exam” statement
- **Default preset was incorrect, operator did not notice**
 - Gain, depth, other settings incorrect
 - Still sufficient to rule out full bladder but did not meet ‘pass criteria’
- **Incorrect preset persisted into Kidney exam**
 - Intervention: requested that operator check the “verify settings” page

AMOS – Kidney Results (N=2 operators/subjects)

	NUMBER OF INSTANCES (CINE-LOOPS)	PROCEDURE EXECUTION (0-3)	INCLUSION OF TARGET ORGAN VOLUME (0-3)	IMAGING PLANE ORIENTATION (0-3)	SWEEP SPEED ADEQUACY (0-3)	AVERAGE IMAGE QUALITY (0-3)	MEASURABILITY (0-3)	EFFECTIVE IMAGING CONTENT (0-3)	SUCCESS SCORE (0-3)	SUCCESS (Y/N)
	Right - Demo 1	10	3.0	2.0	3.0	2.6	1.8	3.0	2.8	Y
	Right - Demo 2	6	2.3	2.7	2.5	2.5	2.2	2.3	2.4	Y
	Left - Demo 1	6	2.9	2.0	2.5	2.6	1.7	3.0	2.8	Y
	Left - Demo 2	6	2.4	2.3	2.5	2.5	2.1	2.3	2.4	Y

RIGHT KIDNEY INSTANCE (CINE-LOOP #)	PROCEDURE EXECUTION (0-3)	INCLUSION OF TARGET ORGAN VOLUME (0-3)	IMAGING PLANE ORIENTATION (0-3)	SWEEP SPEED ADEQUACY (0-3)	AVERAGE IMAGE QUALITY (0-3)	MEASURABILITY (0-3)	EFFECTIVE IMAGING CONTENTS (0-3)	SUCCESS SCORE (0-3)	SUCCESS (Y/N)
Instance 1	3	3	3	3	3	3	3	3	Y
Instance 2	3	3	3	3	3	3	3	3	Y
Instance 3	3	3	3	3	3	3	3	3	Y
Instance 4 - Doppler	1	1	2	2	2	0	1	1.5	N
Instance 5 - Doppler	2	2	3	2	2	2	2	2	Y
Instance 6 - Doppler	2	2	2	2	2	2	2	2	Y
AVERAGE SCORES	2.3	2.3	2.7	2.5	2.5	2.2	2.3	2.4	Y











Demo 2 – Left Kidney Results (N=1 operator/subject)

RIGHT KIDNEY INSTANCE (CINE-LOOP #)	PROCEDURE EXECUTION (0-3)	INCLUSION OF TARGET ORGAN VOLUME (0-3)	IMAGING PLANE ORIENTATION (0-3)	SWEEP SPEED ADEQUACY (0-3)	AVERAGE IMAGE QUALITY (0-3)	MEASURABILITY (0-3)	EFFECTIVE IMAGING CONTENTS (0-3)	SUCCESS SCORE (0-3)	SUCCESS (Y/N)
Instance 1	3	3	2	3	3	3	3	3	Y
Instance 2	3	3	2	3	3	3	3	3	Y
Instance 3	3	2	2	3	3	3	3	3	Y
Instance 4 - Doppler	2	3	2	2	2	1	2	2	Y
Instance 5 - Doppler	2	2	2	2	2	1	2	2	Y
Instance 6 - Doppler	2	2	1	2	2	1	1	1.5	N
AVERAGE SCORES	2.3	2.3	2.7	2.5	2.5	2.2	2.3	2.4	Y

- **Operator inadvertently inverted image left-right**
 - Does not affect image quality, correctable
 - Reverses action on screen, likely disorienting
- **Complex interface (soft keys) likely explanation**
 - Inadvertent setting change not noticed
 - Subsequent exams affected





Interface Instruction Comparison


Ultrasound Machine

STEP 3 OF 4


VERIFY ULTRASOUND SETTINGS

The screen must appear as shown here. **Verify that:**

- **Live image is displayed** (system not frozen due to inactivity); **parallel lines on image change when probe face touched with gel**
 - If the system is frozen, push **BROWN2** (Freeze).
- **Verify the name: AMOS Demo _GMT Date**
 - If incorrect, follow [End Exam](#) steps B-D, and [Create New Patient](#).
- **Probe is correctly identified as 4C RS**
 - If incorrect, replace probe with the **4C RS** probe.
- **Preset selected is "AMOS1kb-Abdominal"** (preset will show truncated "AMOS1k").
 - If incorrect, push **PINK2** (Application); move **TRACKBALL** to highlight 4C RS probe; push **SET** (Trackball Area); move **TRACKBALL** to highlight correct preset; push **SET** (Trackball Area).
- **"Yellow V" (Screen Indicator) is on the left side.**
 - If incorrect, push **BROWN2** (Freeze), then **PURPLE3 Down** (Left/Right), then **BROWN2** (Freeze).
- **Depth is in 12-15 cm range** (14 in this example)
 - Push **BLUE2** Up to reduce, **BLUE2** Down to increase
- **Focus marker is at 5-9 cm range** (8.5 in this example)
 - Push **PURPLE5** Right (Move Up) or **Left** (move Down).



PROCEDURE MILESTONE: UPON COMPLETION OF THE SYSTEM IS READY TO BEGIN THE EXAMINATION




Tap to select the Preset.

Scroll to "Abdomen" and tap Select.

READY TO SCAN!



Bladder Module Acceptability Ratings



Kidney Module Acceptability Ratings

confidence

supported success

all info included

videos helped


images helped

clear

details included

info flow

easy to use



strongly disagree disagree neither agree strongly agree

■ demo 1 □ demo 2

Conclusion

- Initial AMOS demonstrations highly successful
 - All objectives met
 - Bladder and kidney images = most acceptable to excellent
 - Crew feedback = excellent
- Procedure weakness exposed – related to modular nature
- Complex interface can negatively affect performance
- AMOS platform promising as autonomy-enabling tool
 - Exploration
 - Ideal for gradual integration into operations / gain confidence
 - Crew, Surgeon, Ops community, Programs

Recommending further AMOS testing and development

- **Near term: Repeated performance of AMOS kidney and bladder modules by crew with varied background and experience (validation)**
- **Longer term: Additional modules, integration with diagnostic functions, data systems, specific hardware, SME editing, analogs**


Questions?

“Expanding the Boundaries of Space Medicine and Technology”


Backup




Now: ISS Operations



Beyond Low Earth Orbit



Training and autonomous guidance tool



CLINICAL OUTCOME OPTIMIZATION OF RESEARCH

COMf



TAKING A GOOD IMA Kidney and Bladder Ultrasound

Tips for good composition:
To move the optic disc down the subject needs to look up.
To move the optic disc right the subject needs to look right.
In a good composition the optic disc is centered.



ILMA Intubation

PROBE HANDLING: KID
Correct **probe** orientation will ens
Image is not reversed on the screen
Position the probe with the **red**
up.
Slide away from the spine horizon
Use right and left kidney.
Use small probe rotations and tilt
the image to the longest section t



STEP 1 OF 9 ILMA^h Ultrasound Guided IV Access




Diagnosis and Treatment

STEP 3 OF 10: ACQUI
Keep the cross section of the v
the midline and the spine hori
long view, optimize the view ac
two parallel lines toward the m
As you survey the vein in the lo
ABOVE in the region selected fo
site.
You will need to remember the
the needle insertion because y
surface of the skin before nee



SIDE NOTE
For actual patient conta
gloves should be worn i
universal precaution.



STEP 3 OF 10: ACQUI
Keep the cross section of the v
the midline and the spine hori
long view, optimize the view ac
two parallel lines toward the m
As you survey the vein in the lo
ABOVE in the region selected fo
site.
You will need to remember the
the needle insertion because y
surface of the skin before nee



SIDE NOTE
For actual patient conta
gloves should be worn i
universal precaution.



EVALUATED DICOM INSTANCE GROUPS	NUMBER OF INSTANCES (CINE-LOOPS)	PROCEDURE EXECUTION	INCLUSION OF TARGET ORGAN	IMAGING PLANE ORIENTATION	SWEEP SPEED ADEQUACY	AVERAGE IMAGE QUALITY	MEASURABILITY	EFFECTIVE IMAGING CONTENT	SUCCESS SCORE	SUCCESS (Y/N)
		(0-3)	VOLUME (0-3)	(0-3)	(0-3)	(0-3)	(0-3)	(0-3)	(0-3)	
		3	3	3	3	3	3	3	3	Y
		3	3	3	3	3	3	3	3	Y
		3	3	3	3	3	3	3	3	Y
TOTAL NUMBERS AND AVERAGE SCORES		12	3	3	3	3	3	3	3	Y

Demo 1 – Right Kidney Results (N=1 operator/subject)

RIGHT KIDNEY INSTANCE (CINE-LOOP #)	PROCEDURE EXECUTION (0-3)	INCLUSION OF TARGET ORGAN VOLUME (0-3)	IMAGING PLANE ORIENTATION (0-3)	SWEEP SPEED ADEQUACY (0-3)	AVERAGE IMAGE QUALITY (0-3)	MEASURABILITY (0-3)	EFFECTIVE IMAGING CONTENTS (0-3)	SUCCESS SCORE (0-3)	SUCCESS (Y/N)
Instance 1 (0001)	3	2	1	3	2	0	2	2	Y
Instance 2 (0002)	3	2	2	3	3	1	2	2	Y
Instance 3 (0003)	3	2	1	3	2	2	2	2	Y
Instance 4 (0004)	3	2	2	3	3	2	3	2.5	Y
Instance 5 (0005)	3	2	2	3	2	1	2	2	Y
Instance 6 (0006)	3	1	2	3	1	0	1	1	N
*Instance 7 (0037)	3	3	2	3	3	3	3	3	Y
*Instance 8 (0038)	3	3	3	3	3	3	3	3	Y
*Instance 9 (0039)	3	3	3	3	3	3	3	3	Y
*Instance 10 (0040)	3	3	2	3	3	3	3	3	Y
AVERAGE SCORES	3	3	2	3	2.6	1.8	3	2.8	Y

Demo 1 – Left Kidney Results (N=1 operator/subject)

TARGET ORGAN; ACTIVITY; INSTANCE (CINE-LOOP #)	PROCEDURE EXECUTION (0-3)	INCLUSION OF TARGET ORGAN VOLUME (0-3)	IMAGING PLANE ORIENTATION (0-3)	SWEEP SPEED ADEQUACY (0-3)	AVERAGE IMAGE QUALITY (0-3)	MEASURABILITY (0-3)	EFFECTIVE IMAGING CONTENTS (0-3)	SUCCESS SCORE (0-3)	SUCCESS (Y/N)
Kidney Left Instance 1 (0010)	3	2	2	2	3	2	2	2	Y
Instance 2 (0011)	3	2	1	1	2	0	1	1.5	N
Instance 3 (0012)	2	1	1	1	2	0	1	1.5	N
Instance 4 (0013)	3	3	3	2	3	3	2	2.5	Y
Instance 5 (0014)	3	2	3	3	3	2	3	2.5	Y
Instance 6 (0015)	3	3	3	3	3	3	3	3	Y
AVERAGE SCORES	2.8	3	2.1	2	2.6	1.6	3	2.8	Y